



implox Pty.
Ltd.
HEALTHCARE

EN Laerdal Compact Suction Unit (LCSU 3)
Directions for Use

REF Cat No. 88 00 50 & 88 00 60

CAUTION- Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

ES Unidad de succión compacta Laerdal (LCSU 3)
Indicaciones de uso

REF Cat. nº 88 00 50 y 88 00 60

PRECAUCIÓN: La legislación federal de los EE. UU. restringe la venta de este dispositivo a médicos o por prescripción médica.

FR Unité d'aspiration compacte Laerdal (LCSU 3)
Mode d'emploi

REF Cat N° 88 00 50 et 88 00 60

ATTENTION - En vertu de la loi fédérale américaine, la vente de cet appareil n'est autorisée que par un médecin ou sur ordonnance de ce dernier.

DE Laerdal Kompaktes Absauggerät (LCSU 3)
Bedienungsanleitung

REF Kat.-Nr. 88 00 50 & 88 00 60

VORSICHT - Nach US-Bundesgesetzen darf dieses Gerät nur von einem Arzt bzw. auf Anordnung eines Arztes verkauft werden.

IT Aspiratore compatto Laerdal (LCSU 3)
Istruzioni d'uso

REF N. cat. 88 00 50 & 88 00 60

ATTENZIONE -La legge federale statunitense limita la vendita di questo dispositivo ai medici o su loro prescrizione.

NL Laerdal Compact Suction Unit (LCSU 3)
Gebruiksaanwijzing

REF Cat.nr. 88 00 50 & 88 00 60

ATTENTIE- De federale wetgeving in de Verenigde Staten schrijft voor dat dit apparaat uitsluitend mag worden verkocht of voorgeschreven door een arts.

PT Aparelho compacto de aspiração Laerdal
(LCSU 3) Instruções de uso

REF N° de catálogo 88 00 50 e 88 00 60

ATENÇÃO As leis federais (EUA) exigem que este aparelho seja vendido por médicos ou mediante prescrição médica.

SV Laerdal Compact Suction Unit (LCSU 3)
Anvisningar

REF Serie 88 00 50 & 88 00 60

VAR FÖRSIKTIG - Enligt federal lag (U.S.A.) får denna anordning endast säljas av läkare eller på läkarordination.

FI Kannettava Laerdal-imulaite (LCSU 3)
Käyttöohjeet

REF Tuotenrot 88 00 50 ja 88 00 60

HUOMAUTUS: Yhdysvaltain liittovaltion lain mukaan tämän laitteen saa myydä ainoastaan lääkäri tai lääkärin määräyksestä.

DA Laerdal kompakt sug (LCSU 3)
Brugsvejledning

REF Katalognr. 88 00 50 & 88 00 60

FORSIGTIG - I henhold til amerikansk lov må denne anordning udelukkende sælges af eller på opfordring af en læge.

NO Laerdal kompakt sugenhet (LCSU 3)
Bruksanvisning

REF Kat. nr. 88 00 50 & 88 00 60

FORSIKTIG- Føderale lover (USA) begrenser dette utstyret til salg fra eller etter ordre fra lege.

PL Mały ssak elektryczny prod. Laerdal Medical
(LCSU 3) Instrukcja użycia

REF Nr kat. 88 00 50 & 88 00 60










PRZESTROGA - Prawo federalne USA ogranicza sprzedaż tego urządzenia do sprzedaży przez lekarzy oraz na ich zamówienie.

EN English	2
ES Español	11
FR Français	20
DE Deutsch	29
IT Italiano	38
NL Nederlands	47
PT Português	56
SV Svenska	65
FI Suomalainen	74
DA Dansk	83
NO Norsk	92
PL Polish	101

TABLE OF CONTENTS

IEC Symbols	2
Important Safeguards	3
Introduction and Intended Use Statement.....	3
Your Suction Unit System.....	4
Accessory/Replacement Items	4
Set-Up	5
Battery Connection.....	5
How to Operate Your Suction Unit.....	6
Battery Charging	7
Battery Replacement.....	7
Cleaning Instructions.....	7
Troubleshooting.....	8
Specifications/Classifications.....	8
Two-Year Limited Warranty	9
Declaration of Conformity.....	10

IEC SYMBOLS

- | | | |
|--|---|---|
|  Attention, consult directions for use |  Alternating current |  Standby - On/Off |
|  Type BF applied part |  Direct current |  Center positive polarity indicator |
|  Single Use | | |
|  IP12: Vertically falling water drops shall have no harmful effects when the enclosure is tilted at an angle up to 15° on either side of the vertical. | | |
|  The device contains electrical and/or electronic equipment that must be recycled per EC Directive 2002/96/EC - Waste Electrical and Electronic Equipment (WEEE) | | |

IMPORTANT SAFEGUARDS/INTRODUCTION

IMPORTANT SAFEGUARDS

When using electrical products, especially when children are present, basic safety precautions should always be followed. Read all instructions before using. Important information is highlighted by these terms:

DANGER – Urgent safety information for hazards that will cause serious injury or death.

WARNING – Important safety information for hazards that might cause serious injury.

CAUTION – Information for preventing damage to the product.

NOTE – Information to which you should pay special attention.

READ ALL INSTRUCTIONS BEFORE USING.

SAVE THESE INSTRUCTIONS

DANGER

To reduce the risk of electrocution:

1. Do not place or store product where it can fall or be pulled into a tub or sink and come in contact with water.
2. Do not place in or drop into water or other liquid.
3. Do not reach for a product that has fallen into water. Unplug the mains supply cable immediately.

WARNING

To reduce the risk of burns, electrocution, fire or injury to persons:

1. Close supervision is necessary when this product is used by, on, or near children or physically challenged.
2. Use this product only for its intended use as described in this guide.
3. Never operate this product if:
 - a. It has a damaged power cord or plug.
 - b. It is not working properly.
 - c. It has been dropped or damaged.
 - d. It has been dropped into water.
 Return the product to an authorized Laerdal Medical Service Center for evaluation and repair.
4. Keep the power cord away from heated surfaces.

DANGER

This suction unit is a vacuum suction device designed for the collection of nonflammable fluid materials in medical applications only. Improper use during medical applications can cause injury or death. For all medical applications:

1. All suctioning should be done in strict accordance with appropriate procedures that have been established by a licensed medical authority.
2. Some attachments or accessories may not fit the tubing supplied. All attachments or accessories should be checked prior to use to assure proper fit.

INTRODUCTION

Your suction unit is a compact medical suctioning device designed for reliable, portable operation. This suction unit is ideal for providing emergency suction in the field, transport, and hospital environment. Two collection container options include the 800 ml disposable canister and the 300 ml single-use sealed disposable canister. To maximize product life and performance, follow recommended operating and maintenance procedures.

Intended Use Statement

The device is to be used to remove fluids from the airway or respiratory support system. The device creates a negative pressure (vacuum) that draws fluids through disposable tubing that is connected to a collection canister. The fluids are trapped in the collection canister for proper disposal.

YOUR SUCTION UNIT SYSTEM


Inspect the suction unit and all parts before use.

⚠ WARNING-Do not attempt to use if any parts are damaged or missing.

88 00 50 Model (refer to Figures A1 and A2)

1. Bacteria Filter (non-sterile)
2. Connection Tubing (for use w/800 ml Canister) *
3. Display Panel (top of unit)
4. Unit Connection Elbow (for use w/ 800 ml)
5. Vacuum Inlet Port (side)
6. Tethered Plug
7. Wire Canister Bracket (for use w/800 ml canister only)
8. 800 ml Disposable Canister with Lid
9. 90° Canister Connection Elbow
10. Vacuum Regulator Knob (on side)
11. DC Power Input (on side)
12. Battery Door
13. Unit Carry Handle/Catheter Holder
14. 1,8 m (6') Patient Tubing (not shown)
15. AC to DC Adapter/Charger and Line Cord (not shown)
16. High Capacity Rechargeable Battery (not shown) 12V DC Ni-MH
17. Carry Bag with Shoulder Strap (not shown)


*Use only Laerdal replacement part 88 49 05



Figures A1 and A2 - 88 00 50 model

88 00 60 Model (refer to Figures B1 and B2)

1. Display Panel (top of unit)
2. Vacuum Inlet Port (side)
3. 300 ml Disposable Canister with 0.9m (3') Patient Tubing and Internal Bacteria Filter/Fluid Shut-off
4. Canister Tubing Fitting for Disposal
5. Tethered Plug
6. 12V DC Power Input (on side)
7. Vacuum Regulator Knob
8. Battery Door
9. Unit Carry Handle/Catheter Holder
10. AC to DC Adapter/Charger and Line Cord (not shown)
11. High Capacity Rechargeable Battery (not shown), 12V DC Ni-MH



Figures B1 and B2 - 88 00 60 model

ACCESSORY/REPLACEMENT ITEMS

The following items can be purchased separately as accessories or replacements for your suction unit.

Both Models	
Description	Part No.
External Battery Charger (not shown)	88 00 70 04
High Capacity Rechargeable Battery, 12V DC Ni-MH	88 00 70 05
Tethered Plug	88 00 70 06
Replacement Battery Door	88 00 70 07
AC to DC Adapter/Charger	88 44 00
Line Cord USA	88 44 01
Hospital Grade Line Cord USA	88 44 02
Line Cord UK	88 44 03
Line Cord EU	88 44 04
12V DC Power Cord	88 45 00
Vacuum Gauge Test Fixture (not shown)	88 50 00

Model 88 00 50	
Description	Part No.
Unit Connection Elbow	88 00 50 02
800 ml Disposable Canister Pack (incl Canister, Lid, Filter, Canister Elbow, Patient Tubing & Unit Elbow)	88 00 50 04
Wire Canister Bracket (for 800 ml Canister)	88 00 50 06
Carry Bag - w/Shoulder Strap 800 ml	88 46 00
800 ml Disposable Canisters (48 each)	88 47 01
800 ml Disposable Canisters (6 each)	88 47 03
Bacteria Filter (12 each for the 800 ml Canister)	88 49 01
1,8 m (6') Patient Tubing (for 800 ml Canister)	88 49 03
90° Canister Connection Elbow (for 800 ml Canister)	88 49 04
Connection Tubing (6 each for 800 ml Canister)	88 49 05
Model 88 00 60	
Description	Part No.
Carry Bag - w/Strap 300 ml	88 00 60 05
300 ml Disposable Canister with Patient Tubing (10 each)	88 60 00

SET-UP

Battery Connection

NOTE-All models of LCSU 3 are shipped with the battery in place but not connected. Follow the instructions below:

1. Using a coin or straight-blade screwdriver rotate latch to unlocked position (FIG. 1).
2. Remove door by pulling up on latch (FIG. 2).
3. Remove battery from compartment and plug connector into circuit board (FIG. 3).
4. Replace battery and door; rotate latch to locked position.
5. Fully charge battery for 5 hours before using (please see Battery Charging).



Figure 1



Figure 2



Figure 3

WARNING

Do not attempt to connect any type of suction tubing directly into the vacuum inlet port. Only use with a Laerdal approved canister.

NOTE-Always have a spare 300 ml canister within reach in case the canister in use is full or the filter becomes wet.

88 00 50 Model (800 ml Configuration):

1. Place canister into wire canister bracket; ensure inlet port marked <Patient> is accessible.
2. Attach one end of the connection tubing to connection elbow (FIG. 4).
3. Insert unit connection elbow into side port (FIG. 5).
4. Attach other end of connection tubing to side of filter marked <Out>. **NOTE**-Use only Laerdal supplied replacement filters (FIG. 6).
5. The 90° canister elbow connects to the clear side of filter marked <In> and to the top of canister lid marked <Vacuum> (FIG 6). **NOTE**-Verify clear side of filter marked <In> faces canister.
6. Connect the 1.8m (6') patient tubing to canister lid at inlet port labeled <Patient>.
7. Ensure all connections are secure to prevent leakage in the canister/tubing system.
8. Occlude the suction tube and set suction level according to local protocol before suctioning the patient (FIG 7). Be aware that it might be necessary to adjust the suction level during use.




Figure 4




Figure 5



Figure 6



Figure 7

88 00 60 Model (300 ml Configuration):

1. This single-use sealed disposable canister includes an internal filter and fluid shut-off that automatically stops suctioning when it becomes wet (FIG. 8).
2. Firmly attach canister by pushing connection fitting straight into the open side port of unit. (Both the top connection fitting and the bottom support tab secure the canister) (FIG. 9).
3. Securely attach appropriate suction tip to the tubing.
4. Occlude the suction tube and set suction level according to local protocol before suctioning the patient (FIG 7). Be aware that it might be necessary to adjust the suction level during use.

NOTE-While suction unit will continue to operate if tipped on its side, canister capacity will be reduced. Keep a spare 300 ml disposable canister within reach.



Figure 8



Figure 9

HOW TO OPERATE YOUR SUCTION UNIT

Control Panel Symbols (Figure 10)






-  On/Off
-  External power: Supplied from AC to DC Adapter/Charger or 12V DC Power Cord. Illuminates in GREEN when power is supplied.
-  Battery charging: Illuminates in YELLOW. The light will go out when the battery is fully charged.
-  Low battery: Illuminates in RED when battery reaches a discharged state.
-  Suction level setting: Scale/strength illuminates in GREEN. This scale shows the level of suction strength in mmHg.



Figure 10

Power Source Options

AC OPERATION - Plug the 90 degree power connector of AC to DC Adapter/Charger into DC power input (FIG 11) and connect the line cord. Plug the other end of AC to DC Adapter/Charger into a grounded AC supply.

12V DC OPERATION - Plug the small 90 degree power connector of 12V DC power cord into DC power input (FIG 11). Plug large end of cord into 12V DC power receptacle of vehicle.

BATTERY OPERATION - Your unit is equipped with a high capacity rechargeable battery. For initial charge on new unit, fully charge the battery for a minimum of 5 hours (please see Battery Charging).

NOTE—To operate unit from the rechargeable battery, ensure that no external power sources are plugged into the suction unit.

NOTE—During charging or operating, the power supply may become warm to touch; this is normal.



Figure 11

WARNING

If you get the Low Battery Warning symbol, immediately switch to an external power source to avoid an interrupted suction procedure. If the unit does not receive external power or the battery does not get recharged immediately, the low battery indicator light will remain on and the performance of the unit will drop off rapidly and then shut down.

How To Adjust The Vacuum Level

1. Once power source is selected, turn the unit on by pressing the “On” button. The GREEN light, representing external power, will remain lit when external power is connected.
2. Occlude (block) the patient end of the tubing, then adjust vacuum level from 50 to 550 mmHg by turning the vacuum regulator knob clockwise to increase and counter-clockwise to decrease the vacuum (FIG 12). Release and occlude once more to confirm setting. The desired level of vacuum can be viewed on the LED display (FIG 10).

NOTE—The LEDs have two brightness levels. As the vacuum level is adjusted the LEDs will illuminate in progression. When an LED is at half brightness, it indicates that the vacuum level is halfway between the previous fully lit LED and the half brightness LED. *EXAMPLE: If the 150 mmHg LED is fully illuminated and the 200 mmHg LED is at half brightness, this indicates that the suction level is 175 mm Hg. When the 200 mmHg LED illuminates at full brightness, this indicates the unit has reached 200 mmHg. Attention should be taken when setting the different vacuum levels; and it might be necessary to adjust the suction level to local protocol during use.*

3. Connect suction tip or catheter as appropriate.

NOTE—If the unit does not maintain vacuum, refer to Troubleshooting.

NOTE—For 88 00 50 Model, suction ceases when liquid level reaches float shut-off valve located on underside of 800 ml canister lid.
For 88 00 60 Model, suction ceases when liquid level reaches filter located inside 300 ml canister.

NOTE—Dispose of canister and/or contents according to local protocol.

CAUTION—Further suctioning attempts with a full canister may cause damage to the vacuum pump and voids warranty. Equipment service is required if fluid content is aspirated back into the unit.



Figure 12

BATTERY CHARGING

Models 88 00 50 & 88 00 60 Series are equipped with a factory-installed high capacity rechargeable battery. Located on the display panel is the low battery and charge indicator light (FIG 10).

1. Connect the unit to either an AC or DC power source.
2. The green external power light shall be illuminated. The yellow charge indicator will remain lit while the battery is charging.
3. Ensure that the yellow charging light is illuminated when charging begins. As the battery nears a full charge, the yellow battery charging light may flash for several minutes. This is normal. If your unit does not hold a charge, check that the yellow light turns on when external power is applied with the power button "Off". If problems persist, contact an authorized Laerdal Medical Service Center.

NOTE—Recharging the battery to full capacity may take up to 5 hours depending on the depth of discharge.

If unit is not in use for extended periods, the battery should be recharged every 3-6 months. A fully charged battery will provide approximately 45-60 minutes of continuous operation at zero vacuum level (free flow). Unit can be left on charge when not in use.

CAUTION—Completely discharging the battery will shorten the battery life. Do not operate the unit for more than a few minutes if the low battery indicator light is lit. Recharge the battery as soon as possible.

BATTERY REPLACEMENT (refer to Battery Connection in Set Up)

1. Using a coin or straight-blade screwdriver rotate latch to unlocked position.
2. Remove door by pulling up on latch.
3. Remove battery from compartment and unplug connector from circuit board.
4. Install new battery by reversing the above steps.
5. Dispose of battery properly by following local protocol.

CLEANING INSTRUCTIONS

Collection Canister:

1. To remove canister, push power button to turn unit off. Wait for vacuum level to drop.
2. Disconnect external power source from input receptacle on unit (if applicable).
3. Remove canister from unit or holder by disconnecting the elbow, tubing and filter as needed.

NOTE—Insert tethered plug into side port of unit.

4. The 800 ml disposable collection canister and lid are for single-patient use only and must be discarded after use. The 300 ml disposable collection canister is meant for single-use only and must be discarded after use.

NOTE—Before disposal of the 300ml canister, attach free end of tubing to the fitting at the bottom of the canister (FIG 13). This prevents liquid from leaking out of the canister.



Figure 13

Suction Unit:

1. With the power "Off," disconnect the unit from all external power sources.
2. Wipe the outside housing with a clean damp cloth and detergent.

CAUTION—Do not submerge suction unit in water as this will result in damage to vacuum pump.

3. If disinfection is desired, follow the disinfectant manufacturer's recommended instructions and dilution rates carefully.

Tubing:

1. Disconnect tubing and discard; both patient tubing and connection tubing are considered single-patient use only.

Carry Bag:

1. Wipe the bag with a clean, damp cloth soaked with a mild detergent.
2. If disinfection is desired, follow the disinfectant manufacturer's recommended instructions and dilution rates carefully.

Changing Filter (800 ml disposable canister):

1. Filter can be used for up to 2 months of use, but change filter immediately if contamination or discoloration is observed.
2. Remove the bacteria filter by disconnecting it from the suction unit and lid assembly.
3. Replace it with a new Bacteria Filter # 88 49 01 (12 each) and remount it to the suction unit and disposable canister lid.

NOTE—Verify clear side of filter marked <In> faces canister (FIG 14).

NOTE—Use only the Bacteria Filter provided by Laerdal Medical or one of its Distributors. Substitution may lead to contamination of the unit and/or poor performance and will void warranty.



Figure 14

TROUBLESHOOTING

NOTE—Before you return unit to an authorized Laerdal Medical Service Center, follow Troubleshooting protocol below:

⚠ DANGER

Electric shock hazard. Do not attempt to open or remove cabinet, there are no user-serviceable internal components. If service is required, return the suction unit to an authorized Laerdal Medical Service Center. Opening or tampering with the unit will void the warranty.

Problem	Action
Unit does not power on. (Green external power indicator should be illuminated when power is applied.)	<ol style="list-style-type: none"> 1. Check power sources and connections. 2. Ensure wall outlet is live by plugging in a lamp. 3. If operating from 12V DC, ensure DC outlet is live by plugging in known working device such as a cell phone charger. 4. Verify that battery is properly installed, connected and fully charged before use.
Pump runs, but no vacuum.	<ol style="list-style-type: none"> 1. Check that all tubing is connected properly. 2. Check tubing connections for breaks, leaks, or occlusions. 3. Ensure that 800 ml suction canister float shut-off is not activated or that 300 ml canister filter is not clogged. 4. Check for leaks or cracks in canister assembly.
Low vacuum.	<ol style="list-style-type: none"> 1. Use vacuum adjustment knob to increase vacuum level (return to local protocol level after test). 2. Check system for leaks.
Battery will not hold a charge. (Charge indicator should be illuminated if wbattery is connected during charge mode.)	<ol style="list-style-type: none"> 1. Verify that charge light turns on. 2. Check electrical connections during charging. 3. Ensure wall outlet is live by plugging in a lamp.
Battery seems insufficient, does not hold charge.	<ol style="list-style-type: none"> 1. Perform the following test to determine if battery replacement is necessary: <ol style="list-style-type: none"> a. Charge battery as directed. b. Disconnect charging accessory and operate the LCSU 3 at free flow (no suction load and tubing unobstructed) for 20 minutes. <p>If LCSU 3 stops before completing the 20 minutes, contact an authorized Laerdal Medical Service Center for advice regarding battery replacement.</p>

NOTE—If problem is not resolved, contact your authorized Laerdal Medical Service Center:

SPECIFICATIONS/CLASSIFICATIONS

Size - H x W x D inches (cm)

88 00 50 Model (with 800ml collection canister and holder)..... 8.5" x 7.75" x 9" (21.6cm x 19.7cm x22.9cm)

88 00 60 Model (with 300ml collection canister) 7.1" x 10.5" x2.9" (18cm x 26.7cm x 7.4cm)

Weight - lb. (kg)

88 00 50 Model 3.75 lbs. (1.70 kg)

88 00 60 Model 3.375 lbs. (1.53 kg)

Electrical Requirements..... 100-240 VAC 47-63 Hz 0.75 A max; 12 VDC, 33 W max

International Travel - The suction unit is equipped with an AC to DC Adapter/Charger allowing operation on any AC voltage (100-240 VAC, 50/60 Hz). However the correct power cord must be used to connect to adaptable wall power. (See Accessory List on page 4.) **NOTE**-Check power cord for adaptability before using.

Internal Rechargeable Battery 12 VDC

Vacuum Range

88 00 50 & 88 00 60 Model50 to 550 mm Hg (+/- 27.5 mm Hg)

Air Flow @ pump inlet: 27 LPM (free flow) typical (may be less when running from internal battery)

Collection Canister Capacity

88 00 50 Model Disposable 800 ml (cc) Maximum

88 00 60 Model Single-Use Sealed Disposable 300 ml (cc) Maximum

NOTE—If either unit is operated off vertical or on an uneven surface, the collection canister overflow shutoff may activate prematurely, shutting off suction before canister achieves full capacity. Always have a spare replacement canister readily available.

Environmental Conditions

Operating Temperature Range.....	32°F (0°C) - 104°F (40°C)
Operating Relative Humidity.....	0-95%
Operating Atmospheric Pressure.....	10.2 Psi (70 kPa) - 15.4 Psi (106 kPa)
Storage & Transport Temperature Range.....	-40°F (-40°C) - 158°F (70°C)
Storage & Transport Relative Humidity.....	0-95%
Storage & Transport Atmospheric Pressure.....	7.3 Psi (50 kPa) - 15.4 Psi (106 kPa)

Limited Warranty

88 00 50 & 88 00 60 Model.....	Two-years limited, excluding internal battery and collection canister
Internal Battery	90 Days

Certifications

IEC 601-1; IEC 68; CAN/CSA-C22.2 No. 601.1-M90; UL 2601-1, CE EN 60601-1-2, ISO10079-1:1999

Meets RTCA/DO-160E (for battery operation only; commercial aircraft, airborne equipment) DO-160E - section 21 Category M

NOTE - This unit complies with electromagnetic compatibility standards as defined in the included Declaration of Conformity.

Equipment Classifications

With respect to protection from electric shock	Class I and internally powered
Degree of protection against electric shock	Type BF Applied Parts
Degree of protection against ingress of liquids.....	IP12 and standard power supply
Mode of Operation	Intermittent Operation: 30 minutes on, 30 minutes off

WARNING-Do not use equipment in the presence of a flammable or anesthetic gas mixture.

ISO Classification

88 00 50 & 88 00 60 Models - Electrically powered medical suction equipment for field and transport use according to ISO 10079-1:1999 High Flow/High Vacuum

TWO-YEAR LIMITED WARRANTY

The compressor portion of the Laerdal Compact Suction Unit 88 00 50 & 88 00 60 Models (excluding internal rechargeable battery and collection canister) is warranted to be free from defective workmanship and materials for a period of two years from date of purchase. Internal rechargeable batteries are warranted for 90 days. Refer to the Laerdal Global Warranty for terms and conditions on www.laerdal.com. Any defective part(s) will be repaired or replaced at Laerdal Medical's option if the unit has not been tampered with or used improperly during that period. Make certain that any malfunction is not due to inadequate cleaning or failure to follow the instructions. If repair is necessary, contact your authorized Laerdal Medical Service Center for instructions.

NOTE - Be sure to retain a dated proof of purchase document to verify unit is within 2-year warranty period.

NOTE - This warranty does not cover providing a loaner unit or compensating for costs incurred in rental while said unit is under repair.

THERE IS NO OTHER EXPRESS WARRANTY. IMPLIED WARRANTIES, INCLUDING THOSE OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE, ARE LIMITED TO THE DURATION OF THE EXPRESS LIMITED WARRANTY AND TO THE EXTENT PERMITTED BY LAW ANY AND ALL IMPLIED WARRANTIES ARE EXCLUDED. THIS IS THE EXCLUSIVE REMEDY AND LIABILITY FOR CONSEQUENTIAL AND INCIDENTAL DAMAGES UNDER ANY AND ALL WARRANTIES ARE EXCLUDED TO THE EXTENT EXCLUSION IS PERMITTED BY LAW. SOME STATES DO NOT ALLOW LIMITATIONS ON HOW LONG AN IMPLIED WARRANTY LASTS, OR THE LIMITATION OR EXCLUSION OF CONSEQUENTIAL OR INCIDENTAL DAMAGES, SO THE ABOVE LIMITATION OR EXCLUSION MAY NOT APPLY TO YOU.

This warranty gives you specific legal rights, and you may also have other rights which vary from state to state.



DECLARATION OF CONFORMITY

A NOTE FROM LAERDAL

Thank you for choosing a Laerdal Compact Suction Unit. We want you to be a satisfied customer. If you have any questions or comments, please send them to our address on the back cover.

For service, call your authorized Laerdal Medical Service Center:

Phone _____

Purchase Date _____

Serial # _____

DECLARATION OF CONFORMITY

Manufacturer: Sunrise Medical
Address: dba DeVilbiss Healthcare
100 DeVilbiss Drive
Somerset, PA 15501-2125 USA

Product Designation: Laerdal Compact Suction Unit 3 (LCSU 3)

Type/Model: 88 00 50 & 88 00 60 Models

We herewith declare that the above-mentioned product complies with the requirements of EC Directive 93/42/EEC and the following:

Class: IIa, Rule 2

Quality System Standards Applied: ISO 13485:2003

Notified Body: TÜV NORD

MDD: Annex II Applied

Safety Standards Applied: UL 2601-1
EN 60601-1
CAN/CSA 22.2 No. 601.1-M90
ISO 10079-1:1999

EMC Compliance to: IEC 801-2 through 5
CISPR 11/Level B
EN 60601-1-2
RTCA DO-160E

For Battery Operation Only: RTCA/DO-160E,
Section 21 Category M
Commercial Aircraft, Airborne Equipment

Authorized Representative: Sunrise Medical Ltd.
Sunrise Business Park
High Street
Wollaston, West Midlands DY8 4PS
ENGLAND
44-138-444-6688



Units 22-24
60-66 Richmond Road,
Keswick 5035, South Australia
Telephone: 8351 1455
Facsimile: 8293 7377
Email: CustomerService@implox.com
www.implox.com