

**Surgical Hand Immobilizers with anti-tilt Tab**  
**Lead-free and Latex-free**

**Alumi-Hand is sterile and ready to use**

With the affordable Alumi-Hand you trade the hazards of lead for a sterile, disposable, X-ray compatible hand immobilizer. The Alumi-Hand is made with a pliable aluminium interior for easy hand positioning and a foam exterior for patient comfort. Available in large and extra-large sizes and now in a paediatric size.



**Dura-Hand – Easy to sterilize and is reusable**

Just like the Alumi-Hand, the Dura-Hand is lead-free and X-ray compatible. The Dura-Hand is made of pliable aluminium for easy positioning. Available in large and extra-large sizes.

| Product Code | Description   | Qty per Box |
|--------------|---|-------------|
| ISI-759      | Pedi Alumi-Hands 18cm x 17cm<br>Sterile Single Use        | 4           |
| ISI-760      | Large Alumi-Hands 24cm x 21cm<br>Sterile Single Use       | 4           |
| ISI-761      | Extra-Large Alumi-Hands 33cm x 26cm<br>Sterile Single Use | 3           |
| ISI-DH-100   | Pedi DuraHands 18cm x 17cm<br>Reusable                    | 1           |
| ISI-DH-225   | Large DuraHands 24cm x 21cm<br>Reusable                   | 1           |
| ISI-DH-350   | Extra-Large DuraHands 33cm x 26cm<br>Reusable             | 1           |