

# Sakamoto 'Kiyoko' Bathing Care Manikin

**implox** Pty.  
Ltd.  
HEALTHCARE



The lower back smoothly bends anteriorly and posteriorly. The model can be seated in a wheelchair for simulation with a feeling similar to seating a real patient.

The head can be turned laterally to 30 in both directions as well as upward and downward so that the model can take nearly the same postures as the human body.

Since the lower limbs can be freely moved back and forward or abducted, the legs can be spread sufficiently. Practice of assistance in excretion using a bedpan or a urinal can be performed with a true to life feeling.

Various joints including the shoulders, elbows, wrists, hips, knees and ankles move naturally as those of the human body and enable smooth execution of practice of changing clothes and bedding.



- Can be immersed in water.
- Water spontaneously flows into the model so that it can be dipped into water. It drains quickly when lifted from water.
- Tolerates the use of showers, sponge and detergents.

Height: 160cm  
Weight: 13kg  
Material: Vinyl chloride

Part Number	Description
IMPM100-2	Sakamoto 'Kiyoko' Bathing Care Manikin

**implox** Pty.  
Ltd.  
HEALTHCARE

**Units 22-24**  
**60-66 Richmond Road,**  
**Keswick 5035, South Australia**  
**Telephone: 8351 1455**  
**Facsimile: 8293 7377**  
**Email: [CustomerService@implox.com](mailto:CustomerService@implox.com)**  
**[www.implox.com.au](http://www.implox.com.au)**