



- The limbs have seamless structures.
- The head moves in every direction. (Please hold the model by supporting the head as you would do in holding a neonate.)
- The weight of the head has the same proportion as that of a neonate.
- The anterior fontanelle and ribs can be palpated.
- Water infused into the mouth flows out from the urethra so that changing diapers can be practiced.
- Realistic newborn manikin with weak neck muscles

1. Bathing



2. Baby exercise and feeding



3. Dressing / undressing and changing diapers



4. Anal thermometry and various measurements



Height: 49cm
Weight: 3kg
Material: Silicone

Part Number
IMPM107-3
IMPM107-4

Description
Sakamoto Baby Male Version
Sakamoto Baby Female Version

implox Pty.
Ltd.
HEALTHCARE

Units 22-24
60-66 Richmond Road,
Keswick 5035, South Australia
Telephone: 8351 1455
Facsimile: 8293 7377
Email: CustomerService@implox.com
www.implox.com.au