


Sakamoto Transparent Suction Model

implox Pty.
Ltd.
HEALTHCARE



Facilitation of training in Oral, Nasal, Pharyngeal, and Tracheal Suction Techniques.

Relevant Anatomic Structures are Visible.

1. Internal Structure Visible

Observation of tongue, epiglottis, trachea and esophagus. Facilitates understanding of internal anatomy.



2. Verification of tube location by using catheter with LED light.

Catheter with LED light supplied with model for observation of catheter placement.



Catheter Position at Pyriform Sinus



bronchial catheter position

Note: Not Appropriate for Use with Liquids

Includes:

- 1 x Suction Model
- 1 x Suction Tube with LED Light
- 1 x Lubricant Oil
- 1 x Storage Case

Weight: 5.5kg
Case size/W68×D34×H32cm

Sakamoto Transparent Suction Model

implox Pty.
Ltd.
HEALTHCARE

1. Flexion / Extension of Head and Neck



2. Lateral Movement of the Head Part to Left and Right



3. Confirmation of Tube Placement



Part Number
IMPM175-1

Description
Sakamoto Transparent Suction Model

implox Pty.
Ltd.
HEALTHCARE

Units 22-24
60-66 Richmond Road,
Keswick 5035, South Australia
Telephone: 8351 1455
Facsimile: 8293 7377
Email: CustomerService@implox.com
www.implox.com.au