

Sakamoto Suction Simulator

implox Pty.
Ltd.
HEALTHCARE




This model is a simulator with a lifelike texture.

Various suction methods can be practised.

The suction of secretions not only in the oral/nasal cavities but also in tubes inserted via the oral cavities is possible.

Since the trachea has already been incised, suction with tracheal cannula insertion can also be performed.

The head can be moved to the left or right



Includes:

- 1 x Suction Simulator
- 1 x Baby Powder
- 1 x Lubricant Oil
- 1 x Container for Making Artificial Sputum
- 1 x Agitator Stick for Mixing Artificial Sputum
- 1 x Syringe with Tube
- 1 x Silicone Cap
- 1 x Flushing Funnel
- 1 x Flushing Tube
- 1 x Tray

Height: W68cm x D34cm x H32cm
Weight: 15kg

Part Number	Description
IMPM175	Sakamoto Suction Simulator
IMPM	Artificial Sputum

implox Pty.
Ltd.
HEALTHCARE

Units 22-24
60-66 Richmond Road,
Keswick 5035, South Australia
Telephone: 8351 1455
Facsimile: 8293 7377
Email: CustomerService@implox.com
www.implox.com.au

Various Types of Suction Training are Possible



Suction of secretions in the oral cavity



Suction of secretions in inserted tubes



Suction of secretions in the nasal cavity



Suction of secretions in the tracheal cannula



Washing of the bronchi



Washing of the oral cavity