

Sakamoto Baby Touch Vital Signs Simulator

implox Pty.
Ltd.
HEALTHCARE



A High Quality Simulator for Accurate Assessment of Vital Signs in the Neonate



Easy to Read Display



Simple to Use Touch Panel

Set Vital Signs (Heart Rate, Pulse Rate, Respiration Rate, Body Temperature) at multiple levels

Includes:

- 1 x Sakamoto Baby Touch Manikin
- 1 x Simulator
- 1 x Set of Cables
- 1 x Baby Cloth
- 1 x Stethoscope
- 1 x Thermometer
- 1 x Baby Powder

Sakamoto Baby Touch Vital Signs Simulator

implox Pty.
Ltd.
HEALTHCARE

Heart Rate



Auscultate Heart Sounds



Configurable to 6 levels:
(60, 80, 120, 140, 160, 200) bpm

Pulse Rate

Palpation of Pulse at 6 Points, configurable to 6 levels: (60, 80, 120, 140, 160, 200)



Superficial Temporal Arteries



Carotid Arteries



Radial Arteries

Respiration Rate

Realistic Up and Down Movement of Abdominal Area



Configurable to 6 Levels: (20, 40, 50, 60, 80, 100)

Body Temperature



Configurable to 5 levels: 32C - 42C

Part Number	Description
IMPM179	Sakamoto Baby Touch Vital Signs Simulator

Weight/7.2kg (Including accessories)
Case size/W55×D32×H37cm

implox Pty.
Ltd.
HEALTHCARE

Units 22-24

**60-66 Richmond Road,
Keswick 5035, South Australia**

Telephone: 8351 1455

Facsimile: 8293 7377

Email: CustomerService@implox.com

www.implox.com.au