

Sakamoto Musclemate 2 Intramuscular Injection Simulator

implox Pty.
Ltd.
HEALTHCARE



Correct or incorrect injection site is indicated with built-in lamps and a buzzer.



Correct

The green lamp is lit when the injection site and the depth are correct.



Incorrect

The red light is flashed, and the buzzer sounds when the injection site is incorrect.



Incorrect Depth

Both the red and green lamps are flashed and the buzzer sounds when the injection site is correct but too deep.

*The light of the LED lamps is visible through the skin.

An intramuscular injection and palpation simulator.

- Realistic rendition of human skin and underlying muscle.
- Indicates the correct or incorrect injection site including the depth with LED lamps.

Preparation



Palpation



Pinching



Deltoid Muscle



Subcutaneous Injection



implox Pty.
Ltd.
HEALTHCARE

Units 22-24

**60-66 Richmond Road,
Keswick 5035, South Australia**

Telephone: 8351 1455

Facsimile: 8293 7377

Email: CustomerService@implox.com

www.implox.com.au

Part Number
IMPM118-1

Description
Sakamoto Musclemate Intramuscular Injection Simulator