

# Sakamoto Intramuscular Injection Model Demonstration Type II

**implox** Pty.  
Ltd.  
HEALTHCARE



The buzzer and flashing red light warn in case of incorrect points.

Shibusanbu Method, Clarke's and Hochstetter's Point Identification are Facilitated in the Model.

- Inbuilt skeletal structure facilitates identification of correct injection site by palpation of the iliac crest and greater trochanter.
- Correct needle insertion is indicated by illumination of a green light. the red light illuminates and a buzzer is activated when the incorrect injection site is used or if insertion is too deep. Fluid can be injected.



Injection training to correct injection sites is possible.



Fluid can be injected.

**implox** Pty.  
Ltd.  
HEALTHCARE

**Units 22-24**

**60-66 Richmond Road,  
Keswick 5035, South Australia**

**Telephone: 8351 1455**

**Facsimile: 8293 7377**

**Email: [CustomerService@implox.com](mailto:CustomerService@implox.com)**

**[www.implox.com.au](http://www.implox.com.au)**

Part Number  
IMPM152-1

Description  
Sakamoto IM Injection Model Demonstration Type 2