


Sakamoto Gluteal Injection Trainer

implox Pty.
Ltd.
HEALTHCARE



A lateral recumbent model for injections at Clarke's point and Hochstetter's point.

- Being a lateral recumbent position, the model permits more practical training.
- The epidermis is made of a unique material that leaves no needle mark.



Hochstetter's Point



Clarke's Point

The skeletal landmarks such as the anterior superior iliac spine, posterior superior iliac spine, and greater trochanter can be felt, and the model is ideal for practicing the determination of the injection site.

Part Number
IMPM152-3

Description
Sakamoto Gluteal Injection Trainer

implox Pty.
Ltd.
HEALTHCARE

Units 22-24
60-66 Richmond Road,
Keswick 5035, South Australia
Telephone: 8351 1455
Facsimile: 8293 7377
Email: CustomerService@implox.com
www.implox.com.au