


# Sakamoto IM Injection Model Demonstration

**implox** Pty.  
Ltd.  
HEALTHCARE



- Inbuilt skeletal structure facilitates identification of correct injection site by palpation of the iliac crest and greater trochanter.
- Correct needle insertion is indicated by illumination of a green light. the red light illuminates and a buzzer is activated when the incorrect injection site is used or if insertion is too deep. Fluid can be injected.



Part Number  
IMPM152

Description  
Sakamoto IM Injection Model Demonstration

**implox** Pty.  
Ltd.  
HEALTHCARE

**Units 22-24**  
**60-66 Richmond Road,**  
**Keswick 5035, South Australia**  
**Telephone: 8351 1455**  
**Facsimile: 8293 7377**  
**Email: [CustomerService@implox.com](mailto:CustomerService@implox.com)**  
**[www.implox.com.au](http://www.implox.com.au)**