


Sakamoto Buttockmate IM Injection Simulator with Alarms

implox Pty.
Ltd.
HEALTHCARE



The epidermis is made of a unique material that leaves no needle mark.



Selection of an incorrect site is indicated by blinking of the red LED and buzzing sound. Successful puncturing at the correct site is indicated by the green lamp. (Only at Clarke's Point and Hochstetter's Point)



Being compact, the device needs a small space for storage.

Reproducing Clarke's Point and Hochstetter's Point.

Simulation for Intramuscular Injection the Lateral Recumbent Position.



Skeletal reference points such as the anterior superior and anterior posterior iliac spines and greater trochanter can be palpated for realistic training in the identification of the correct injection site.

Fluid can be injected.

implox Pty.
Ltd.
HEALTHCARE

Units 22-24

**60-66 Richmond Road,
Keswick 5035, South Australia**

Telephone: 8351 1455

Facsimile: 8293 7377

Email: CustomerService@implox.com

www.implox.com.au

Part Number
IMPM153-1

Description
Sakamoto Buttockmate Type 2 IM Injection Simulator