

Sakamoto IV Trainer with Indicator



A “Fit-on” Simulator for enhanced training in IV access technique.

Built in lamps, audible feedback and a buzzer indicate correct or incorrect injection into the three main veins: Basilic vein, Cephalic vein and Median cubital vein


Built-in lamps, chime sound and a buzzer indicate correct or incorrect positioning.



Correct Placement
Green lamp, chime sound



Incorrect Positioning
Red lamp, buzzer



Blood circulated by an electric pump (included)
Observe flashback

Weight: 4kg
Case: W51×D42×H20cm

Sakamoto IV Trainer with Indicator

implox Pty.
Ltd.
HEALTHCARE



Intravenous Injection



Blood Sampling



Insertion of Indwelling Needle / Catheter



Fixing with Adhesive Tape



Part Number	Description
IMPM183	Sakamoto Trainer IV Arm with Indicator
IMPM182	Arm for Simulator (M183)
IMPM183-1	Exchangable Blood Vessels Pk of 100

implox Pty.
Ltd.
HEALTHCARE

Units 22-24
60-66 Richmond Road,
Keswick 5035, South Australia
Telephone: 8351 1455
Facsimile: 8293 7377
Email: CustomerService@implox.com
www.implox.com.au