



implox

IMPLOX PTY. LTD.
 Medical Equipment Suppliers
 Units 22-24
 60-68 Richmond Road
 Keswick 5035
 South Australia

Telephone: (08) 8351 1455
 Facsimile: (08) 8293 7377
 email: sales@implox.com
 Web: www.implox.com

2018 NEW OFFERING ENVENOMATION MANAGEMENT BANDAGE

APRIL, 2018

The AEROFORM Snakebite Bandage

A Bite Injury inflicted by venomous reptiles requires rapid emergency intervention with a reliable and easily used compression bandage.

The Australian Resuscitation Council has reported that poor outcomes in the management of snakebite injury may be related to inadequate Pressure Immobilisation Technique (ARC, 9.4.8, 2011)

The innovative AEROFORM Snakebite Bandage carries an indicating strip made up of a series of rectangles. When correct pressure is applied during bandaging, the rectangles "**reform**" into squares, telling the rescuer that **correct pressure** has been applied.

The AEROFORM Snakebite Bandage is an essential addition to any First Aid Kit!

Available in standard (adult) size (**10.5 metre stretched**) or small (child) size (**4.5 metre stretched**)

Prices do not include freight and delivery



AFHSB100 (12 x Standard Bandages)

10.5m x 10cm

AFHSB100/1 (1 x Standard Bandage)

10.5m x 10cm

AFSB10 (12 x Small Bandages)

4.5m x 10cm

AFBS10/1 (1 x Small Bandage)

4.5m x 10cm