



B.I.G Bone Injection Gun

Leading intraosseous device



1. Four different locations to inject give more options depending on the situation.
2. Not relying on external power (no batteries involved) - the device is as ready after 5 years as it was on the first year. Other devices need to be monitored to know what the level of battery is.
3. No need for capital purchase for each caregiver/ambulance - the needle of the B.I.G and the device are the one and the same unit.
4. Sterility: the cover of the B.I.G prevent the needle from being contaminated before injected into the bone. (An exposed needle cannot guarantee sterility).
5. Compact size - small and effective.
6. No need to clean device, single use to be discarded after use saving time and money.
7. Linear penetration - no micro fractures and body reaction. The way it penetrates prevents mistakes in case of movements by the patient or care taker.
8. The BIG has superior safety elements for the care taker- the protection of the needles by the casing prevents any accidental puncture to the paramedic/ doctor.
9. Stability - the B.I.G. allows proper stability of the needle, no matter what the penetration depth is.
10. Price - the B.I.G. device supplied complete, sterile and ready to use making it extremely price competitive.