

seca

Precision for health

Flat Scales



seca added value

Flat scales for heavyweights with a separate display element? There is such a thing – from seca of course. It is the platform and bariatric scales seca 634. The advantage of the mobile display is obvious. Medical staff can read the results from any location – discreetly and conveniently at a desk, at eye-level on the wall or from a stand in the middle of the room.



It's easy to step onto the extra large, slip-proof platform.



The patient's BMI can be determined when his height is entered.

seca 634

Platform and bariatric scales

Obesity is the fastest growing health problem in the world. This battery-operated digital scale has been especially developed for weighing obese patients, with a weighing range of up to 360 kg and a non-slip 56 x 56 cm platform, only 5.5 cm high – a sturdy base for overweight or obese patients. A chair may also be used on the scale and will not influence the weighing result due to the pre-TARE function. The cable remote display (cable length: 2.5 m) with easy-to-read LCD can be held, wall mounted or placed on a table. The 100 g graduation means that the smallest changes in weight are registered, even after the patient has left the scale – thanks to the auto-HOLD function. The integrated Body Mass Index (BMI) function also emphasises the many uses of this unique multifunctional scale.

- Capacity: 360 kg / 800 lbs
- Graduation: 100 g / 0.2 lbs
- Functions: TARE, pre-TARE, auto-HOLD, BMI, kg/lbs switch-over, damping, automatic switch-off
- Optional: stand for cable remote display seca 472, RS232 interface seca 460, power adapter seca 447, switch-mode power adapter seca 400



Accessories

seca 472

Stand for cable remote displays of
seca scales

The seca 472 stand is the stable solution for cable remote displays of seca scales seca 634, seca 656 and seca 674. The sturdy construction with non-slip steel feet can be positioned everywhere quickly and securely. The display element can easily be attached to the stand by means of a stable bracket.

seca added value

A mobile weighing and measuring system? There is such a thing – from seca of course. With the addition of the adapter element seca 437 for the measuring rod seca 217, the flat scales seca 874 turn into a mobile and very stable 2-in-1-solution (seca 254, page 36). Weighing and measuring can be done at the same time – with the same quality as with permanently installed equipment.

The integrated handle and low weight make these scales easy to transport.



The angle of the display can be adjusted for easy reading and can be used separately from the base.

seca 869

Flat scales with cable remote display

A two-meter-long cable connects the display with the weighing base. Advantage: the display can be set up separately from the base and read from its position on the wall or desk. Additional convenience is supplied by the many extra weighing functions: BMI for nutritional conditions, HOLD to keep measurements displayed and TARE for weighing small children held in an adult's arms. It is suitable for stationary or mobile use, thanks to its handle, compact size, low weight and a leveling base for set-up on all types of surfaces.

- Capacity: 250 kg / 550 lbs
- Graduation: 100 g < 150 kg > 200 g / 0.2 lbs < 330 lbs > 0.5 lbs
- Functions: TARE, HOLD, BMI, kg/lbs switch-over, automatic switch-off
- Optional: stadiometer seca 217 combined with adapter element seca 437, power adapter seca 447, switch-mode power adapter seca 400, carrying case seca 414 (for combination with stadiometer and adapter element) and seca 421



The two push buttons on the front side of the scale can be operated with the foot.

seca 874

Mobile flat scales for mobile use with push buttons and double display

The seca 874 can do everything that a scale in mobile use should do and more. It features low weight, extremely robust construction and easy-to-use design with long-life batteries and precision unaffected by humidity levels. When conducting one examination after the other, there's no need to bend over repeatedly, thanks to the push button that can be turned on and off with a toe tap. The double display lets the patient and medical personnel read the results from two different prospects at the same time. And thanks to the 2 in 1 button, the scale can weigh a child held in its parent's arms.

- Capacity: 200 kg / 440 lbs
- Graduation: 50 g < 150 kg > 100 g / 0.1 lbs < 330 lbs > 0.2 lbs
- Functions: mother/child-function, auto-HOLD, kg/lbs switch-over, automatic switch-off
- Optional: stadiometer seca 217 combined with adapter element seca 437, carrying case seca 414 (for combination with stadiometer and adapter element) and seca 421



Both doctor and patient can see the result at the same time in the double display.

The low weight of only 4.2 kilograms makes mobile use possible.



Tip-on function for weighing without turning on the scales.

seca 876

Flat scales for mobile use

Scales intended for mobile use have to be lightweight. The compact seca 876 weighs only 4.2 kilograms. This scale is nevertheless very stable and suitable for persons weighing up to 250 kilograms. It stands firmly on four leveling bases that guarantee stability, even when small children are weighed while held by an adult. The mother/child-function ascertains the tare of the weight with just the touch of a button.

- Capacity: 250 kg / 550 lbs
- Graduation: 100 g < 150 kg > 200 g / 0.2 lbs < 330 lbs > 0.5 lbs
- Functions: mother/child-function, kg/lbs switch-over, tip-on automatic switch-on, automatic switch-off
- Optional: stadiometer seca 217 combined with adapter element seca 437, carrying case seca 414 (for combination with stadiometer and adapter element) and seca 421



The extra-wide surface, 400 x 45 x 300 mm, and the flat shape facilitate trouble-free weighing of heavyweight patients.

seca 813

Electronic flat scales with very high capacity

Because increasing numbers of people weigh more, there are these scales that easily hold up to 200 kilograms on their extra wide platform. The tough rubber coating takes the heaviest challenges lightly and the flat construction makes stepping on very easy.

- Capacity: 200 kg / 440 lbs
- Graduation: 100 g / 0.2 lbs
- Functions: auto-HOLD, step-off, kg/lbs switch-over, tip-on automatic switch-on, automatic switch-off



Even the basic model is complex in design and robustly ribbed for a secure stance during the weighing process.

seca 803

Electronic flat scales with high-quality two-component rubber surface

The classic, clean line of the basic model turns heads and the high-quality rubber coating feels good to the touch. It's also pleasantly easy to read the weights in the display because the LCD figures are unusually large.

- Capacity: 150 kg / 330 lbs
- Graduation: 100 g / 0.2 lbs
- Functions: auto-HOLD, kg/lbs switch-over, tip-on automatic switch-on, automatic switch-off
- Optional: carrying case seca 421

Due to 28-millimeter LCD figures, the weights in the display are easy to read.



Accessories

seca 421

Carrying case



With this carrying case of water-repellent material, the following seca scales can be safely transported: seca 869, seca 874, seca 876 and seca 803.



The face of the dial is especially large to make reading the figures much easier.



The rugged metal housing withstands years of use under the hardest conditions imaginable.

seca 762

Mechanical personal scales with fine 500 g graduation

The precise design classic has a large clear and well-defined dial with a graduation of 500 grams. With a large, low-level, non-slip platform and scratch resistant finish it is easy for the patient to step on and weigh.

- Capacity: 150 kg, 320 lbs, 150 kg / 320 lbs
- Graduation: 500 g, 1 lbs, 500 g / 1 lbs
- Functions: reset-to-zero function



The large dial makes it easy to read results.

seca 750

Robust mechanical floor scales

Pretty and stable: with their pure, light design, these scales will make a good impression in any clinic. A stable, powder-coated steel casing, a circular display and an easy to clean tread area of imitation leather – what more does one need for comfortable weighing? The seca 750 is available in three versions: kilograms, pounds and kilograms/pounds.

- Capacity: 150 kg, 320 lbs, 150 kg / 320 lbs
- Graduation: 1 kg, 1 lbs, 1 kg / 1 lbs
- Functions: reset-to-zero function



The housing is of metal treated with powdercoated lacquer.



The stable material on the platform ensures safety from slippage.



**Units 22-24
60-66 Richmond Road,
Keswick 5035, South Australia
Telephone: 8351 1455
Facsimile: 8293 7377
Email: CustomerService@implox.com
www.implox.com**