



Laryngoscopes Product Catalogue



Index

Green System Fibre Optic Laryngoscope Handles.....	Page 1
Greenspec2 Reusable Fibre Optic Laryngoscope Blades.....	Page 2
GreenLite Disposable Fibre Optic Laryngoscope Blades.....	Page 2
Conventional Illumination Laryngoscope Handles.....	Page 3
Reusable Conventional Illumination Laryngoscope Blades.....	Page 4
Disposable Conventional Illumination Laryngoscope Blades.....	Page 5
Disposable Handle & Blade Combinations.....	Page 6
Laryngoscope Sets.....	Pages 7-11
Difficult Intubation.....	Page 12
Difficult Intubation Training.....	Page 13

Green System Fibre Optic Laryngoscope Handles



GreenSpex

- Perfectly made to ISO 7376(3) internationally recognized specifications
- Fits all "green-system" blades made to this standard
- Xenon bulb gives incomparable cold bright light in excess of 3,000 Lux through fibre-optic blade element
- One-second screw removal of bulb housing allows easy autoclaving
- Smooth barrel knurling provides perfect grip without risk of tearing anesthetist's gloves
- See wide selection of blades + set combinations

4211X	Green Spex Medium Handle
4211	Green Spex Mini Handle
4211X-R	Green Spex Recharge Medium Handle
4212X-R	Green Spex Recharge Mini Handle
4213X	Green Spex Stubby Handle
5072	Trucharge Battery AA Cell



5036AU	Trucharge Kit - Charger Base, Mini + Medium GreenSpex Handles & Batteries
5035AU	Trucharge Kit - Charger Base, 2 x Medium GreenSpex Handles and Batteries
5020AU	Trucharge + Power Supply + Plug & 2 x Medium Handles
5030AU	Trucharge + Power Supply + Plug + Medium & Mini Handle
552501	Trugrip Handle Cover



TruLED Rechargeable Laryngoscope Handle

- Indestructibly engineered from stainless steel 316 that prevents normal wear & tear caused by repeated mounting of blades
- Guaranteed maintenance free for 5 years
- Resistant to rust and chemicals
- Can be autoclaved and sterilized repeatedly by all prevalent systems
- Unique one ideal size ergonomic design with silicon cushioned grip
- Led light source that provides bright bluish white light spectrum
- Unique patented configuration emits uniform and focused light
- Fits all green system blade made to ISO 7376(3) international standards
- Rechargeable Lithium ion battery with latest technology, no capacity deterioration; power for over 10 hours continuous use between charges ; recharge time equals use time; RoHS compliant for banned materials
- Wide range of blade and set combinations available

4220AU	Tru-LED Handle Set, Cradle and Power Supply
--------	---

GreenSpec2 Reusable F.O. Laryngoscope Blades



GreenSpec 2


The most extensive range of blade sizes with Xenon illumination

- Widest choice of Macintosh and Miller Blades and the unique Oxyport paediatric blades
- Incomparable cold, bright light - 2800 LUX
- Simple, trouble-free use. No screws, no lost parts
- Autoclavable FibreClip™ light guide, tested to withstand over 500 cycles
- Stainless steel parts guaranteed for 5 years

- 4240 Green Spec 2 Miller Oxy 00
- 4241 Green Spec 2 Miller Oxy 0
- 4249 Green Spec 2 Miller 00
- 4250 Green Spec 2 Miller 0
- 4251 Green Spec 2 Miller 1
- 4252 Green Spec 2 Miller 2
- 4253 Green Spec 2 Miller 3
- 4254 Green Spec 2 Miller 4

- 4230 Green Spec 2 Macintosh 0
- 4231 Green Spec 2 Macintosh 1
- 4232 Green Spec 2 Macintosh 2
- 4233 Green Spec 2 Macintosh 3
- 4234 Green Spec 2 Macintosh 4
- 4235 Green Spec 2 Macintosh 5

GreenLite Disposable F.O. Laryngoscope Blades



GreenLite Disposable Blade


Uncompromising Performance in the fight against cross infection

- Unique patented stainless steel design exactly duplicates metal blade performance and feel
- Excellent light transmission using highest acrylic light guides
- Provides visualization without “dazzle” effect
- Simple, reliable “click” locking and unlocking mechanism
- Fits onto Green system handles
- Individually packaged colour coded for size and ready to use
- Unique design allows for perfect disposability

- 4319 Green Lite Miller 00
- 4320 Green Lite Miller 0
- 4321 Green Lite Miller 1
- 4322 Green Lite Miller 2
- 4323 Green Lite Miller 3

- 4341 Green Lite Macintosh 1
- 4342 Green Lite Macintosh 2
- 4343 Green Lite Macintosh 3
- 4344 Green Lite Macintosh 4

Conventional Illumination Laryngoscope Handles



Dolphin Handle

Patented time-saving concept that makes cleaning easy with new stainless steel head to prolong life.

- Perfectly made to ISO 7376 : 2009 internationally recognized specifications
- All Stainless Steel parts , Head, barrel and cap
- Fits all conventional blades made to this standard
- Allows immersion in cleaning solutions without removal of batteries
- Mildly knurled barrel provides perfect grip without risk of tearing anesthetist's gloves
- Competitively priced for highest standard of quality
- Withstands repeated autoclaving with removal of batteries
- Available individually and in sets with wide range of blade sizes

5572 Water Proof Mini Handle Rough
5571 Water Proof Standard Rough Handle



UltraSafe Handle

Stainless steel heads providing long-lasting service for conventional handles

- Perfectly made to ISO 7376(1) internationally recognized specifications
- Fits all conventional blades made to this standard
- Stainless steel contact points increase life
- Smooth barrel knurling provides perfect grip without risk of tearing anesthetist's gloves
- Competitive price uncompromised by poor quality
- Wide range of blade and set combinations available

5523 Stubby Rough Handle
5526 Ultra Safe Mini Handle
5527 Ultra Safe Paediatric Handle
5525 Ultra Safe Standard Handle

Reusable Conventional Laryngoscope Blades




Equip Blades

Gold standard for value, Gold standard for quality

- Stainless Steel Medical Grade
- Guaranteed for 5 years
- Each blade has easy to read "manufacturing date"
- Streamlined profile
- Conforms to all ISO standards
- Wide variety of blade sizes and set combinations

5810	Equip Macintosh 0
5811	Equip Macintosh 1
5812	Equip Macintosh 2
5813	Equip Macintosh 3
5814	Equip Macintosh 4
5815	Equip Macintosh 5

5827	Equip Miller Oxy 00
5828	Equip Miller Oxy 0
5829	Equip Miller 00
5830	Equip Miller 0
5831	Equip Miller 1
5832	Equip Miller 2
5833	Equip Miller 3
5834	Equip Miller 4



Shucman 2 F.O. Blades For Use with Conventional Handles

The patented upgrade to long term fibre-optic illumination

- Only fibre-optic blade that fits onto a conventional handle
- Allows upgrade to fibre-optic, xenon bulb illumination (2,400 Lux)
- Patented concept
- Unique replaceable fibre-clip extends blade life indefinitely
- Large variety of blade sizes and set combinations available

5240	Shucman 2 Macintosh 0
5241	Shucman 2 Macintosh 1
5242	Shucman 2 Macintosh 2
5243	Shucman 2 Macintosh 3
5244	Shucman 2 Macintosh 4
5255	Shucman 2 Macintosh 5

5248	Shucman 2 Oxy 00
5249	Shucman 2 Oxy 0
5247	Shucman 2 Miller 00
5250	Shucman 2 Miller 0
5251	Shucman 2 Miller 1
5252	Shucman 2 Miller 2
5253	Shucman 2 Miller 3
5254	Shucman 2 Miller 4

Disposable Conventional Laryngoscope Blades

EquipLite Blades

Disposable stainless steel laryngoscope blade with LED light source. Compatible with conventional system laryngoscope handles.



- Prevents cross-infection
- Saves hidden reprocessing costs
- LED light technology means strong, focused, reliable light
- Ideal blue-white spectrum
- Optimal viewing
- Longer battery life
- Atraumatic tip for patient safety and comfort
- Specially designed for safer, cleaner disposability
- Color coded for easy size identification
- Same great blade profiles in full range of Macintosh and Miller sizes
- Unique design allows for perfect disposability

4819	Equip Lite Miller 00
4820	Equip Lite Miller 0
4821	Equip Lite Miller 1
4822	Equip Lite Miller 2
4823	Equip Lite Miller 3
4824	Equip Lite Miller 4

4840	Equip Lite Macintosh 0
4841	Equip Lite Macintosh 1
4842	Equip Lite Macintosh 2
4843	Equip Lite Macintosh 3
4844	Equip Lite Macintosh 4
4845	Equip Lite Macintosh 5

Lite-Blade

“... only single-use plastic blade that fits conventional handle”



- Ergonomic slim design with tooth protection area
- Feels and performs exactly like a metal blade
- Disposable and “single-patient-use” blades decrease the high risk of cross-infection
- Tested for tensile strength well in excess of the physical forces required during intubation
- Unbeatable, strong, brilliant white light
- Individually packaged colour coded for size and ready to use
- Available in all Macintosh and Miller sizes

4041	Lite Blade Slim Macintosh 1
4042	Lite Blade Slim Macintosh 2
4043	Lite Blade Slim Macintosh 3
4044	Lite Blade Slim Macintosh 4

4020	Lite Blade Miller 0
4021	Lite Blade Miller 1
4022	Lite Blade Miller 2
4023	Lite Blade Miller 3

Disposable Laryngoscope Blade & Handle Combinations



TruLite

- The most preferred disposable stainless steel blade/re-enforced plastic handle combination
- LED illumination gives natural white colour with spectrum now recognized as ideal for enhancing viewing of tissues
- Soldered trip ensures non-traumatic blade entry
- Special light weight, only 120 gram including the AA batteries
- Unbreakable handle that folds so that blade-handle fits into most EMS pouches
- Ideal for Hospital Crash Trolleys
- Macintosh 2,3 & 4 colour coded sizes already available
- Miller 1,2 & 3 sizes available
- Requires 2 x AA Batteries

4642	TruLite Macintosh 2 - One Piece Disposable Laryngoscope
4643	TruLite Macintosh 3 - One Piece Disposable Laryngoscope
4644	TruLite Macintosh 4 - One Piece Disposable Laryngoscope
4619	TruLite Miller 00 - One Piece Disposable Laryngoscope
4620	TruLite Miller 0 - One Piece Disposable Laryngoscope
4621	TruLite Miller 1 - One Piece Disposable Laryngoscope
4622	TruLite Miller 2 - One Piece Disposable Laryngoscope
4623	TruLite Miller 3 - One Piece Disposable Laryngoscope
4624	TruLite Miller 4 - One Piece Disposable Laryngoscope

Laryngoscope Sets



4306 Green Lite EMS Set: Macintosh 2,3,4 Blades & Green Spex Handle



4260/TL-AU Green Spec 2 Set: Tru-Led Handle + Macintosh 1,2,3
4260/X Green Spec 2 Set: Standard X-Handle + Macintosh 1,2,3



4261/TL-AU Green Spec 2 Set: Tru-Led Handle + Macintosh 1,2,3,4
4261/X Green Spec 2 Set: X-Handle + Macintosh 1,2,3,4



4262/X Green Spec 2 Set: X-Handle + Macintosh 2,3,4

Laryngoscope Sets



4470X Green Spec 2 Set: Operating Room
Includes: 'Pro' Blade Macintosh 3.5, Stubby Handle,
Green speX Mini & Medium Handles, Macintosh 1,2,3,4
Blades and Miller Oxy 0, Miller 00 & 1 Blades



5710E Equip Set: Ultrasafe Medium Handle + Macintosh 1,2,3



5711E Equip Set: Ultrasafe Medium Handle + Macintosh 1,2,3,4



5713E Equip Set: Ultrasafe Medium Handle + Macintosh 2,3,4

Laryngoscope Sets



5430D

Shucman 2 Set: Mini Waterproof Handle & Miller 00,0,1



5439D

Shucman 2 Set: Waterproof Handle & Macintosh 2,3,4



5440D

Shucman 2 Set: Waterproof Handle & Macintosh 1,2,3



5441D

Shucman 2 Set: Waterproof Handle & Macintosh 1,2,3,4

Laryngoscope Sets



5570D Shucman 2 Set: Operating Room
Includes 'Pro' Blade Macintosh 3.5, Stubby Handle, Mini Handle, Dolphin Medium Handle, Macintosh 1, 2, 3, 4 & Miller Oxy 0, Miller 00 & 1 Blades



4770X Green-Spec Pro Set 'Pro' 3.5 Macintosh & Greenspex Handle & Spare Xenon Bulb



5443M Tru-MR™ 6 Blade Set - Contains Handle, 6 Blades (Macintosh 2,3,4, Miller 1,2,3) And 2 Battery Packs



4774X Green Spec Pro Set 'Pro' Blade Macintosh 4.5 with Greenspex Handle & Spare Xenon Bulb

Laryngoscope Sets



5770

Shucman-Pro Set 'Pro' Macintosh 3.5 Blade & Dolphin Medium Handle & Spare Bulb



5771

Shucman-Pro Set 'Pro' Macintosh 3.5 Blade & Dolphin Medium Handle & Spare Bulb



5442M

Tru-MR™ 5 Blade Set - Contains Handle, 5 Blades (Macintosh 2,3,4, Miller 0,1) And 2 Battery Packs



- 5161E2S Truview EVO2 Set Shucman 2: 1 x Small Blade, 1 x Mini Handle, 1 x Optishape Stylet
- 4161E2S Truview EVO2 Set: Greenspec 2: 1 x Small Blade, 1 x Mini Handle, 1 x Optishape Stylet

- 5162E2S Truview Maxi EVO2 Set Shucman 2: 1 x Adult Blade, 1 x Small Blade, 1 x Medium handle, 1 x Mini Handle, 2 x OptiShape Stylets
- 4162E2S Truview Maxi EVO2 Set Greenspec 2: 1 x Adult Blade, 1 x Small Blade, 1 x Medium handle, 1 x Mini Handle, 2 x OptiShape Stylets

- 4164E2S Truview EVO2 GreenSpec 2 Set: 1 x Adult Blade, 1 x Medium Handle, 1 x Optishape stylet
- 5164E2S Truview EVO2 Shucman 2 Set: 1 x Adult Blade, 1 x Medium Handle, 1 x Optishape stylet

- 4160A/TL Truview EVO2 GreenSpec 2 Set: 1 x Adult Blade, 1 x Small Blade, 1 x Truled Handle
- 4160E2S Truview EVO2 GreenSpec 2 Set: 1 x Adult Blade, 1 x Small Blade, 1 x Medium Handle
- 5160E2S Truview EVO2 Shucman 2 Set: 1 x Adult Blade, 1 x Small Blade, 1 x Waterproof Handle

- 4163E2S Truview EVO2 GreenSpec Set Infant: 1 x Infant Blade, 1 x Truled Infant Handle, 1 x Infant Optishape stylet
- 4270/MX Truview GreenSpec 2 Set: 1 x Mini Handle, 1 x Miller 00, 0, 1 Blades

- 410001 ETT Guide-Optishape 5-6.5mm
- 410004 ETT Guide-Optishape 7-9mm

Difficult Intubation Training



- 4136AU Truview Premier Intubation Trainer Kit Greenspec 2:
1 x Adult Blade, 1 x Small Blade, 1 x medium Handle, 1 x mini
Handle, 2 x Optishape stylets, 1 x Premier Display, 1 x Carry case
- 5136AU Truview Premier Intubation Trainer Kit Shucman 2:
1 x Adult Blade, 1 x Small Blade, 1 x medium Handle, 1 x mini
Handle, 2 x Optishape stylets, 1 x Premier Display, 1 x Carry case

implox Pty.
Ltd.
HEALTHCARE

**Units 22-24
60-66 Richmond Road,
Keswick 5035, South Australia
Telephone: 8351 1455
Facsimile: 8293 7377
Email: CustomerService@implox.com
www.implox.com**