



For Fast & Safe Evacuations

Safe...

The AlbacMat has been designed with safety, for both the rescuer and patient, being of paramount importance

Effective ...

The ability of the AlbacMat to travel easily over most surfaces including stairs, carpet & concrete – allows for improved efficiency when evacuating a facility

Single person operation...

Evacuations can be performed by one person per patient in most situations which is essential when time is critical

Use in any location...

The AlbacMat is totally independent of its environment. It doesn't matter where that patient may be – in bed, the bathroom or the corridor, they can still be placed onto the AlbacMat and taken to safety

Ready for instant use...

Its compact, self storage and lightweight design means it can be readily available in varying locations for instant use



Design

- Tapered for patient comfort
- Easily accessible dual straps for securing the patient
- Fluorescent yellow straps for detection in dim lighting
- Foot rest pocket for the comfort and security of the patient
- Simple instructions are clearly visible on the self storage pocket
- Three separate handles for individual rescuer's requirements
- Insulated board provides patient protection from uncomfortable surfaces
- Patient can remain on the AlbacMat until the danger has passed and still be treated if required
- Handles provide excellent control on staircases
- Compact design makes storage easy

Product Range

AlbacMat Rescue Mat – Ideal for aged care / residential facilities and general public use (Product No: ALBA01)

AlbacMat Rescue Mat with Thermal Cover – Larger design for specific use in emergency services & mining areas, but is also suitable where the AlbacMat Rescue Mat is used (Product No: ALBA02)

AlbacMat Holder – Holds 4 AlbacMats (either size) for easy storage & quick access (Product No: ALBA03)