



The smart and small AED

Heartstart FR2 provides emergency first responders with an effective means to treat victims in Sudden Cardiac Arrest (SCA). This lightweight, rugged and easy-to-use Automated External Defibrillator (AED) allows first responders to deliver lifesaving therapy as soon as they arrive on the scene. Its patented SMART™ Biphasic waveform is the most patient-friendly technology available today.

**The Heartstart FR2 is...**

**Small and Light...**

The FR2 weighs in at only 2.1 kg, and measures only 218 x 218 x 66 mm, so it is easy both to transport and store.

**Easy to Use...**

Instant response is critical. Both audible and visual prompts guide the user through the simple steps required.

**Smart...**

Utilises clinically proven SMART™ Biphasic technology. Combined with SMART™ analysis system, the FR2 provides safe and effective therapy for the individual patient.

**Features:**

**Long life battery:** Minimum 300 shocks when new. Shelf life before installation and standby life after installation typically 5 years each.

**Rechargeable training battery:** A cost efficient way to utilise the device for training

**Backlit display:** Improves the quality of ECG and text prompts in conditions of direct sunlight and darkness

**Selftest and status indicator:** Monitors readiness and clearly displays whether or not the device is operational

**Easy and flexible configuration:** Allows extensive configuration possibilities, including treatment protocol in accordance with new ERC or AHA guidelines

**Manual Override:** (Model M3840A only) can be set on or off

P/N: M3840A Defi brillator Heartstart EM Model FR2+

P/N: M3841A Defi brillator Heartstart S Model FR2