

Laerdal AED Trainer 2

Realistic, advanced and affordable AED trainer which simulates the HeartStart FR2 AED in appearance, operation and voice prompts.

The Laerdal AED Trainer 2 offers an affordable alternative for educating the layperson and healthcare provider in the effective use of an automatic external defibrillator. This trainer was also developed to specifically simulate the HeartStart FR2 defibrillator

Used together with AED Little Anne or AED Resusci Anne. (see pages..) the AED Trainer 2 with Laerdal Link technology offers unique realism and response to pad placement.

The trainer has 10 easily selectable Guidelines 2005 compatible scenarios, which can be overridden and adjusted from the optional remote control.

Latest software version and computer software for configuration and custom scenarios are available as free downloads from www.laerdal.com

The trainer uses 6 C-size alkaline batteries, included.



940050 Laerdal AED Trainer 2

Laerdal AED Trainer 2 Includes:

Trainer, Carry Case, 6 Batteries, and Standard Adult Training Pads.



945050 Remote Control

945055 Programming Kit

Contains programming cable (RS-232 serial cable) plus CD with latest software in all languages, Configuration and Scenario Builder software. The software on the CD is also available as free download from www.laerdal.com



07-10900 Standard Adult Training Pads

M3871A Standard Pediatric Training Pads

945090 Laerdal LINK Training Pads

These pads are for use with AED Little Anne or AED Resusci Anne



945030 Carry Case AED Trainer 2