

HeartStart FRx Trainer

A realistic training partner to the HeartStart FRx

The HeartStart FRx Trainer simulates how the HeartStart FRx will operate during a real-life sudden cardiac arrest event.

Clear, confident voice instructions provide easy-to-follow directions, including CPR instructions, identical to the ones given by the real HeartStart FRx AED.

Eight easily selectable Guidelines 2005 compatible scenarios provides versatile and realistic training.

Multiple use training pads gives realistic feedback when used on a manikin with either internal or external manikin adapters. The same training pads can be used with the HeartStart FRx AED. Use Infant/Child key to activate Infant/Child mode.

The trainer is powered by 4 AA-size batteries (not included).



861307 Heartstart FRx Trainer

HeartStart FRx Trainer Includes:

Trainer, Training Pads II, 1 External Manikin Adapter strip, Carry Case and Instruction for Use.

Consumables



989803139271 Heartstart Training Pads II
(containing one set of Training Pads II in training pads case)

989803139291 Replacement Training Pads II
(pair of training pads on disposable liner for use in training pads case)

989803139311 Infant/Child Key

M5089A External Manikin Adapter; 5 pack

M5088A Internal manikin adapters.

Self adhesive patches to be placed on the inside of a Little Anne or Resusci Anne manikin skin.

942510 Carry Case



implox Pty.
Ltd.
HEALTHCARE

Units 22-24
60-66 Richmond Road,
Keswick 5035, South Australia
Telephone: 8351 1455
Facsimile: 8293 7377
Email: CustomerService@implox.com
www.implox.com