

## HeartStart Trainer

*An affordable, safe and effective tool for AED training.*

The HeartStart Trainer simulates how the HeartStart HSI will operate during a real-life sudden cardiac arrest event.

Clear, confident voice instructions provide easy-to-follow directions, including CPR instructions, identical to the ones given by the real HSI AED.

Eight easily selectable Guidelines 2005 compatible scenarios provides versatile and realistic training.

Multiple use adult or child training pads gives realistic feedback when used on a manikin with either internal or external manikin adapters. The same training pads can be used with the HeartStart HSI AED.

The trainer is powered by 4 AA-size batteries (not included).



**M5086A** Heartstart Trainer

**HeartStart Trainer includes:**

Trainer; Adult Training Pads, 1 External Manikin Adapter strip, Carry Case and Instruction for Use.

### Consumables



**M5073A** Adult training pads cartridge

**M5074A** Infant/child training pads cartridge

**M5093A** Replacement adult training pads

**M5094A** Replacement infant/child training pads

**M5089A** External manikin adapters 5 pack

**M5088A** Internal manikin adapters.

Self adhesive patches to be placed on the inside of a Little Anne or Resusci Anne manikin skin.



**942500** Carry Case

**implox** Pty.  
Ltd.  
HEALTHCARE

**Units 22-24**  
**60-66 Richmond Road,**  
**Keswick 5035, South Australia**  
**Telephone: 8351 1455**  
**Facsimile: 8293 7377**  
**Email: CustomerService@implox.com**  
**www.implox.com**