

Nursing Kid



Laerdal Simulation Solutions



Laerdal[®]
helping save lives

Nursing Kid

Nursing Kid realistically represents a six year-old child, designed for skills and scenario-based training of care and management of a variety of pediatric hospital patients.

- **Educationally effective** for training pediatric patient care skills
- Programmable scenarios provide **standardized training** while customizable scenarios and real-time instructor control allows training to be adapted to meet student needs
- **Highly mobile** for use in multiple clinical settings and hospital transport

- Head with anatomical landmarks, trachea and esophagus, along with simulated lungs and stomach, allows the practice of many procedures, including NG, OG, tracheal care and suctioning, and insertion, securing and care of endotracheal tubes
- **Programmable sounds** for the auscultation and diagnosis of normal and abnormal heart, breath and bowel sounds
- **Intravenous therapy and site care** using the multivenous pediatric IV arm

Airway Features:

- Irrigation of the eye and ear (simulated)
- Nasal packing
- Ear, eye, nose and mouth care
- Insertion, securing and care of endotracheal tubes
- Insertion and suctioning of oropharyngeal and nasopharyngeal airways
- Tracheostomy care and tracheal suctioning
- NG tube insertion, care and removal
- Gastric lavage and gavage
- Various oxygen delivery procedures
- Manually generated carotid pulse

Catheterization:

- Interchangeable male and female genitalia with connector valves and urinary and colon reservoirs
- Complete urinary catheterization
- Will retain indwelling or straight catheter
- Enema procedures may be performed using fluid for realistic return

Other Features:

- Full range of motion for realistic patient handling
- Dressing and bandaging techniques possible
- Pediatric transfer techniques possible

Sounds - VitalSim™:

- Lung sounds synchronized with breathing rate, 0-60 BPM
- Normal and abnormal heart sounds
- Individual lung or bilateral sound selection
- Normal and abnormal bowel sounds
- Prerecorded or live feed voice sounds (with optional wireless microphone)

- Articulating IV arm with replaceable skin and infusible vein system allows peripheral intravenous therapy and site care
- Venipuncture possible in the antecubital fossa and dorsum of the hand
- Accessible veins include median, basilic and cephalic
- Bilateral thigh, bilateral deltoid, gluteal and ventrogluteal intramuscular injections

VitalSim Unit required for operation.
Sold separately.

350-05050 Nursing Kid (VitalSim Capable)
350-10050 Nursing Kid (Non-VitalSim Capable)

Nursing Kid Includes:

Full-body Pediatric Manikin, Airway Lubricant, Simulated Blood, Hospital Gown, Directions for Use and Carry Case

implox Pty.
Ltd.
HEALTHCARE

Units 22-24

60-66 Richmond Road,
Keswick 5035, South Australia
Telephone: 8351 1455
Facsimile: 8293 7377

Email: CustomerService@implox.com
www.implox.com